




① Unit: TAAC. Air Conditioning Laboratory Unit

Always included in the supply:

Teaching Technique used


SCADA. EDIBON Computer Control System



② Control Interface Box



③ Data Acquisition Board




④ Software for:
- Computer Control
- Data Acquisition
- Data Management



Computer (not included in the supply)


⑤ Cables and Accessories
⑥ Manuals



www.edibon.com

Products
Products range
Units
9. Thermodynamics & Thermotechnics

PROCESS DIAGRAM AND ELEMENTS ALLOCATION



**OPEN CONTROL
+
MULTICONTROL
+
REAL TIME CONTROL**

5 actuators and
16 sensors controlled
from any computer, and
working simultaneously

ST-1 Temperature sensor of dry bulb, type J.
ST-2 Temperature sensor of wet bulb, type J.
ST-3 Temperature sensor of dry bulb, type J.
ST-4 Temperature sensor of wet bulb, type J.
ST-5 Temperature sensor of dry bulb, type J.
ST-6 Temperature sensor of wet bulb, type J.
ST-7 Temperature sensor of dry bulb, type J.
ST-8 Temperature sensor of wet bulb, type J.
ST-9 Temperature sensor (evaporator outlet), type J.

ST-10 Temperature sensor (condenser outlet), type J.
ST-11 Temperature sensor (evaporator inlet), type J.
ST-12 Temperature sensor (steam generator).
SC-1 Refrigerant flow sensor.
SP-1 Pressure sensor (outlet of the condenser).
SP-2 Pressure sensor (inlet of the condenser).
SP-3 Differential pressure sensor (measure of flow).
DSP-1 Bourdon manometer (outlet of condenser).
DSP-2 Bourdon manometer (inlet of evaporator).

DSP-3 Bourdon Manometer (outlet of evaporator).
AP-1 High-pressure cut-out.
AN-1 Level switch.
AR-1 Steam generator resistance.
AR-2 Pre-heater.
AR-3 Re-heater.
Compressor.
Fan with speed control.

Items supplied as standard

① TAAC. Unit:

This unit has as objective to introduce the student in the world of the air conditioning installations, as well as to study and determine the good parameters for the unit operation in function of the environmental demands (humidity, heat, temperature and refrigeration, etc).

Anodized aluminium structure and panels in painted steel.

Main metallic elements in stainless steel.

Diagram in the front panel with similar distribution that the elements in the real unit.

Tunnel of 300 x 300 x 1600 mm. , made in stainless steel with 2 windows of 200 x 300 mm. to visualize the tunnel inside.

2 Electrical heating resistances (computer controlled): one of 2000W (pre-heater) to the inlet of the evaporator and other of 1000 W (re-heater) to the outlet of the evaporator.

4 Hygrometers placed along the tunnel, formed each one by 2 temperature sensors (wet and dry bulb).

Fan, with speed control from computer, 0.25KW, 2500 r.p.m, Qmax 2160 m³/h.

Evaporator.

Compressor, 1/2 CV.

Condenser unit. At 5°C = 1591W. 980 m³/h.

High-pressure cut-out, tared at 14 bar.

Filter dryer.

Flow meter and refrigerant flow sensor, range: 5-60 l./h.

Temperature sensors (11):

4 dry bulb "J" type, 4 wet bulb "J" type, 1 "J" type (inlet of the evaporator), 1 "J" type (outlet of the evaporator), 1 "J" type(outlet of the condenser). Sensors range: -40 to 750° C.

Pressure sensors (3):

1 sensor (outlet of the condenser), 1 sensor (inlet of the condenser), 1 differential sensor (measure of flow).

Bourdon manometers (3):

1 bourdon manometer (outlet of the condenser), 1 bourdon manometer (inlet of the evaporator), 1 bourdon manometer (outlet of the evaporator).

Psychometric chart and Enthalpy diagram of R134a.

The unit incorporates wheels for mobility.

②TAAC/CIB. Control Interface Box :

Control interface box with process diagram in the front panel and with the same distribution that the different elements located in the unit, for an easy understanding by the student.

All sensors, with their respective signals, are properly manipulated for -10V. to +10V computer output. Sensors connectors in the interface have different pines numbers (from 2 to 16), to avoid connection errors. Single cable between the control interface box and computer.

The unit control elements are permanently computer controlled, without necessity of changes or connections during the whole process test procedure. Simultaneously visualization in the computer of all parameters involved in the process. Calibration of all sensors involved in the process.

Real time curves representation about system responses. Storage of all the process data and results in a file. Graphic representation, in real time, of all the process/system responses.

All the actuators' values can be changed at any time from the keyboard allowing the analysis about curves and responses of the whole process. All the actuators and sensors values and their responses are placed in only one computer screen.

Shield and filtered signals to avoid external interferences.

Real time PID control with flexibility of modifications from the computer keyboard of the PID parameters, at any moment during the process. Real time PID and on/off control for pumps, compressors, resistances, control valves, etc. Real time PID control for parameters involved in the process simultaneously. Proportional control, integral control and derivative control, based on the real PID mathematical formula, by changing the values, at any time, of the three control constants (proportional, integral and derivative constants). Open control allowing modifications, at any time and in a real time , of parameters involved in the process simultaneously.

Possibility of automatization of the actuators involved in the process.

Three safety levels, one mechanical in the unit, other electronic in control interface and the third one in the control software.

③ DAB. Data Acquisition Board:

PCI Data acquisition board (National Instruments) to be placed in a computer slot. Bus PCI.

Analog input:

Number of channels= 16 single-ended or 8 differential. Resolution= 16 bits, 1 in 65536.

Sampling rate up to: 250 KS/s (Kilo samples per second).

Input range (V) = ±10V. Data transfers=DMA, interrupts, programmed I/O. DMA channels=6.

Analog output:

Number of channels=2. Resolution= 16 bits, 1 in 65536. Maximum output rate up to: 833 KS/s.

Output range(V) =±10V. Data transfers=DMA, interrupts, programmed I/O.

Digital Input/Output: Channels=24 inputs/outputs. D0 or DI Sample Clock frequency: 0 to 1 MHz.

Timing: Counter/timers=2. Resolution: Counter/timers: 32 bits.

④TAAC/CCSOF. Computer Control+Data Acquisition+Data Management Software:

Compatible with actual Windows operating systems. Graphic and intuitive simulation of the process in screen. Compatible with the industry standards.

Registration and visualization of all process variables in an automatic and simultaneously way. Flexible, open and multicontrol software, developed with actual windows graphic systems, acting simultaneously on all process parameters. Analog and digital PID control. Menu for PID and set point selection required in the whole work range. Management, processing, comparison and storage of data. Sampling velocity up to 250,000 data per second guaranteed.

Student calibration system for all sensors involved in the process.

It allows the registration of the alarms state and the graphic representation in real time.


Comparative analysis of the obtained data, after to the process and modification of the conditions during the process. Open software, allowing to the teacher to modify texts, instructions. Teacher's and student's passwords to facilitate the teacher's control on the student, and allowing the access at different work levels.

This unit allows that the 30 students of the classroom can visualize simultaneously all results and manipulation of the unit, during the process, by using a projector.

⑤Cables and Accessories, for normal operation.

⑥Manuals: This unit is supplied with 8 manuals: Required Services, Assembly and Installation, Interface and Control Software, Starting-up, Safety, Maintenance, Calibration & Practices Manuals.

***References 1 to 6: TAAC + TAAC/CIB + DAB + TAAC/CCSOF + Cables and Accessories + Manuals are included in the minimum supply, enabling a normal operation.**




TAAC Unit



TAAC/CIB



DAB



TAAC/CCSOF

Complementary items to the standard supply

PLC. Industrial Control using PLC (7 and 8):

⑦ PLC-PI. PLC Module:

Circuit diagram in the front panel.

Front panel:

Digital inputs(X) and Digital outputs (Y) block:

16 Digital inputs, activated by switches and 16 LEDs for confirmation (red).

14 Digital outputs (through SCSi connector) with 14 LEDs for message (green).

Analog inputs block:

16 Analog inputs (-10V. to + 10V.) (through SCSi connector).

Analog outputs block:

4 Analog outputs (-10V. to + 10V) (through SCSi connector).

Touch screen:

High visibility and multiple functions.

Display of a highly visible status.

Recipe function.

Bar graph function.

Flow display function.

Alarm list.

Multi language function.

True type fonts.

Back panel:

Power supply connector.

Fuse 2A.

RS-232 connector to PC.

Inside:

Power supply outputs: 24 Vdc, 12 Vdc, -12 Vdc, 12 Vdc variable.

Panasonic PLC:

High-speed scan of 0.32 μ sec. for a basic instruction.

Program capacity of 32 Ksteps, with a sufficient comment area.

Free input AC voltage(100 to 240 V AC).

DC input: 16 (24 V DC).

Relay output: 14 (250 VA AC/2 A).

High-speed counter.

Multi-point PID control.

Digital inputs/outputs and analog inputs/outputs Panasonic modules.

Communication RS232 wire, to computer (PC).



PLC-PI

⑧ TAAC/PLC-SOF. PLC Control Software:

For this particular unit, always included with PLC supply.


Items available on request

⑨ TAAC/CAL. Computer Aided Learning Software (Results Calculation and Analysis).

⑩ TAAC/FSS. Faults Simulation System.

Software Main Screens

Main screen



Note: ST=Temperature sensor. SP= Pressure sensor. AR=Heating resistance.

Examples of Sensors Calibration screens

Analog Input Channel: AI 8

Least Squares Fit

Gain: 96.615 Offset: 0.626275

Points to average: 0

Volts: 0.111 Calibrated: 9.7651

ENTER DONE

R134a Flow

Min Value: 0

Max Value: 60

ENTER DONE

Simultaneous Calibration

Reference Value: 0 Data taken: 0

Reference Select	Sensors	Volts	Calibrated	ΔT
<input type="checkbox"/>	ST-1	0.2753	29.8346	29.83
<input type="checkbox"/>	ST-2	0.3335	29.7856	29.79
<input type="checkbox"/>	ST-3	0.3291	29.0641	29.06
<input type="checkbox"/>	ST-4	0.3254	29.5453	29.55
<input type="checkbox"/>	ST-5	0.3295	29.4276	29.43
<input type="checkbox"/>	ST-6	0.3458	34.752	34.75
<input type="checkbox"/>	SC-1	0.0037	0.0037	0
<input type="checkbox"/>	SC-2	0.004	0.01826	0.02
<input type="checkbox"/>	SC-3	0.0059	0.0059	0
<input type="checkbox"/>	SC-4	0.0059	0.0059	0
<input type="checkbox"/>	SC-5	0.0059	0.0059	0
<input type="checkbox"/>	SC-6	0.0059	0.0059	0
<input type="checkbox"/>	SC-7	0.0059	0.0059	0
<input type="checkbox"/>	SC-8	0.0059	0.0059	0
<input type="checkbox"/>	SC-9	0.0059	0.0059	0
<input type="checkbox"/>	SC-10	0.0059	0.0059	0
<input type="checkbox"/>	SC-11	0.0059	0.0059	0
<input type="checkbox"/>	SC-12	0.0059	0.0059	0
<input type="checkbox"/>	SC-13	0.0059	0.0059	0
<input type="checkbox"/>	SC-14	0.0059	0.0059	0
<input type="checkbox"/>	SC-15	0.0059	0.0059	0
<input type="checkbox"/>	SC-16	0.0059	0.0059	0
<input type="checkbox"/>	SC-17	0.0059	0.0059	0
<input type="checkbox"/>	SC-18	0.0059	0.0059	0
<input type="checkbox"/>	SC-19	0.0059	0.0059	0
<input type="checkbox"/>	SC-20	0.0059	0.0059	0
<input type="checkbox"/>	SC-21	0.0059	0.0059	0
<input type="checkbox"/>	SC-22	0.0059	0.0059	0
<input type="checkbox"/>	SC-23	0.0059	0.0059	0
<input type="checkbox"/>	SC-24	0.0059	0.0059	0
<input type="checkbox"/>	SC-25	0.0059	0.0059	0
<input type="checkbox"/>	SC-26	0.0059	0.0059	0
<input type="checkbox"/>	SC-27	0.0059	0.0059	0
<input type="checkbox"/>	SC-28	0.0059	0.0059	0
<input type="checkbox"/>	SC-29	0.0059	0.0059	0
<input type="checkbox"/>	SC-30	0.0059	0.0059	0
<input type="checkbox"/>	SC-31	0.0059	0.0059	0
<input type="checkbox"/>	SC-32	0.0059	0.0059	0
<input type="checkbox"/>	SC-33	0.0059	0.0059	0
<input type="checkbox"/>	SC-34	0.0059	0.0059	0
<input type="checkbox"/>	SC-35	0.0059	0.0059	0
<input type="checkbox"/>	SC-36	0.0059	0.0059	0
<input type="checkbox"/>	SC-37	0.0059	0.0059	0
<input type="checkbox"/>	SC-38	0.0059	0.0059	0
<input type="checkbox"/>	SC-39	0.0059	0.0059	0
<input type="checkbox"/>	SC-40	0.0059	0.0059	0
<input type="checkbox"/>	SC-41	0.0059	0.0059	0
<input type="checkbox"/>	SC-42	0.0059	0.0059	0
<input type="checkbox"/>	SC-43	0.0059	0.0059	0
<input type="checkbox"/>	SC-44	0.0059	0.0059	0
<input type="checkbox"/>	SC-45	0.0059	0.0059	0
<input type="checkbox"/>	SC-46	0.0059	0.0059	0
<input type="checkbox"/>	SC-47	0.0059	0.0059	0
<input type="checkbox"/>	SC-48	0.0059	0.0059	0
<input type="checkbox"/>	SC-49	0.0059	0.0059	0
<input type="checkbox"/>	SC-50	0.0059	0.0059	0
<input type="checkbox"/>	SC-51	0.0059	0.0059	0
<input type="checkbox"/>	SC-52	0.0059	0.0059	0
<input type="checkbox"/>	SC-53	0.0059	0.0059	0
<input type="checkbox"/>	SC-54	0.0059	0.0059	0
<input type="checkbox"/>	SC-55	0.0059	0.0059	0
<input type="checkbox"/>	SC-56	0.0059	0.0059	0
<input type="checkbox"/>	SC-57	0.0059	0.0059	0
<input type="checkbox"/>	SC-58	0.0059	0.0059	0
<input type="checkbox"/>	SC-59	0.0059	0.0059	0
<input type="checkbox"/>	SC-60	0.0059	0.0059	0
<input type="checkbox"/>	SC-61	0.0059	0.0059	0
<input type="checkbox"/>	SC-62	0.0059	0.0059	0
<input type="checkbox"/>	SC-63	0.0059	0.0059	0
<input type="checkbox"/>	SC-64	0.0059	0.0059	0
<input type="checkbox"/>	SC-65	0.0059	0.0059	0
<input type="checkbox"/>	SC-66	0.0059	0.0059	0
<input type="checkbox"/>	SC-67	0.0059	0.0059	0
<input type="checkbox"/>	SC-68	0.0059	0.0059	0
<input type="checkbox"/>	SC-69	0.0059	0.0059	0
<input type="checkbox"/>	SC-70	0.0059	0.0059	0
<input type="checkbox"/>	SC-71	0.0059	0.0059	0
<input type="checkbox"/>	SC-72	0.0059	0.0059	0
<input type="checkbox"/>	SC-73	0.0059	0.0059	0
<input type="checkbox"/>	SC-74	0.0059	0.0059	0
<input type="checkbox"/>	SC-75	0.0059	0.0059	0
<input type="checkbox"/>	SC-76	0.0059	0.0059	0
<input type="checkbox"/>	SC-77	0.0059	0.0059	0
<input type="checkbox"/>	SC-78	0.0059	0.0059	0
<input type="checkbox"/>	SC-79	0.0059	0.0059	0
<input type="checkbox"/>	SC-80	0.0059	0.0059	0
<input type="checkbox"/>	SC-81	0.0059	0.0059	0
<input type="checkbox"/>	SC-82	0.0059	0.0059	0
<input type="checkbox"/>	SC-83	0.0059	0.0059	0
<input type="checkbox"/>	SC-84	0.0059	0.0059	0
<input type="checkbox"/>	SC-85	0.0059	0.0059	0
<input type="checkbox"/>	SC-86	0.0059	0.0059	0
<input type="checkbox"/>	SC-87	0.0059	0.0059	0
<input type="checkbox"/>	SC-88	0.0059	0.0059	0
<input type="checkbox"/>	SC-89	0.0059	0.0059	0
<input type="checkbox"/>	SC-90	0.0059	0.0059	0
<input type="checkbox"/>	SC-91	0.0059	0.0059	0
<input type="checkbox"/>	SC-92	0.0059	0.0059	0
<input type="checkbox"/>	SC-93	0.0059	0.0059	0
<input type="checkbox"/>	SC-94	0.0059	0.0059	0
<input type="checkbox"/>	SC-95	0.0059	0.0059	0
<input type="checkbox"/>	SC-96	0.0059	0.0059	0
<input type="checkbox"/>	SC-97	0.0059	0.0059	0
<input type="checkbox"/>	SC-98	0.0059	0.0059	0
<input type="checkbox"/>	SC-99	0.0059	0.0059	0
<input type="checkbox"/>	SC-100	0.0059	0.0059	0

ENTER DONE

Example of PID controls screen

Set Point: 45.00

output range: output high 10.00, output low 0.00

PID gains:

proportional gain (Kc): 1.000


integral time (Ti, min): 0.010

derivative time (Td, min): 0.000

SET

Continue...

Some typical exercises results



○ Actions:

1.-Starting up the steam generator.

● Effects: (changes marked with red dots)

- Temperature is starting to raise in the steam generator.


○ Actions:

1.-Steam generation is going up to 99.5°C.
2.-We start the fan, from position "0" to position "5".

● Effects: (changes marked with red dots)

- Air Flow is raising to 20 litres/second.
- Temperature changes in several ST sensors.
- Pressure changes in SP sensors.
- Flow changes in SC-1 sensor.

-The open control, multicontrol and real time control system is used.




○ Actions:

1.-Steam is generating.
2.-Fan is working now at maximum speed, position "10".
3.-Starting up the compressor.

● Effects: (changes marked with red dots)

- The outlet temperatures in ST-5/ 6/ 7/ 8 sensors are going down.
- The temperatures of the input air in ST-3 and ST-4 sensors are raising.
- Pressure changes in SP sensors.
- Flow changes in SC-1 sensor.

-The open control, multicontrol and real time control system is used.



EXERCISES AND PRACTICAL POSSIBILITIES

Some Practical Possibilities of the Unit:

- 1.- Demonstration of the processes and components used in heating, cooling, humidification, de-humidification of an airstream.
- 2.- Obtaining of the steam generator efficiency curve.
- 3.- Energy balance in the steam generator.
- 4.- Efficiency determination of the preheating resistance.
- 5.- Preheating effect in an air conditioning installation.
- 6.- Dehumidification process study.
- 7.- Material balance in the evaporator.
- 8.- Energy balance in the evaporator.
- 9.- Re-heat effect.
- 10.- Experimental determination of the air specific heating capacity.

Other possible practices:


- 11.- Psychrometric chart.
- 12.- Example of the air properties determination.
- 13.- Usage of psychrometric chart.
- 14.- Determination of the airflow.
- 15.- Temperature sensors calibration.
- 16.- Pressure sensors calibration.
- 17.- Determination of a PWM controller adjustment parameters.
- 18.- Properties of the Refrigerant R134a.
- 19.- Enthalpy-Pressure diagram for the refrigerant R134a.

Practices to be done by PLC Module (PLC-PI) + PLC Control Software:

- 20.- Control of the TAAC unit process through the control interface box without computer.
- 21.- Visualization of all the sensors values used in TAAC unit process.
- 22.- Calibration of all sensors included in TAAC unit process.
- 23.- Hand on of all the actuators involved in the TAAC unit process.

- 24.- Realization of different experiments, in automatic way, without having in front the unit. (This experiment can be decided previously).
- 25.- Simulation of outside actions, in the cases do not exist hardware elements (Example: test of complementary tanks, complementary industrial environment to the process to be studied, etc).
- 26.- PLC hardware general use and manipulation.
- 27.- PLC process application for the TAAC unit.
- 28.- PLC structure.
- 29.- PLC inputs and outputs configuration.
- 30.- PLC configuration possibilities.
- 31.- PLC program languages.
- 32.- PLC different programming standard languages (literal structured, graphic, etc.).
- 33.- New configuration and development of new process.
- 34.- Hand on an established process.
- 35.- To visualize and see the results and to make comparisons with the TAAC unit process.
- 36.- Possibility of creating new process in relation with the TAAC unit.
- 37.- PLC Programming Exercises.
- 38.- Own PLC applications in accordance with teacher and student requirements.

POSSIBILITIES OF OTHER AVAILABLE EXPANSIONS



ORDER INFORMATION

Items supplied as standard

Minimum configuration for normal operation includes:

- ① Unit: TAAC. Air Conditioning Laboratory Unit.
- ② TAAC/CIB. Control Interface Box.
- ③ DAB. Data Acquisition Board.
- ④ TAAC/CCSOF. Computer Control + Data Acquisition + Data Management Software.
- ⑤ Cables and Accessories, for normal operation.
- ⑥ Manuals.

*** IMPORTANT: Under TAAC we always supply all the elements for immediate running as 1, 2, 3, 4, 5 and 6.**

Complementary items to the standard supply

PLC. Industrial Control using PLC. (7 and 8):

- ⑦ PCL-PI. PLC Module.
- ⑧ TAAC/PLC-SOF. PLC Control Software.
- ⑨ TAAC/CAL. Computer Aided Learning Software (Results Calculation and Analysis). (Available on request).
- ⑩ TAAC/FSS. Faults Simulation System. (Available on request).

Expansions

- ⑪ Mini ESN. Multipost EDIBON Mini Scada-Net System.
- ⑫ ESN. Multipost EDIBON Scada-Net System.

REQUIRED SERVICES

Electrical supply: 220V., 1-phase + neutral + ground ,50 Hz. , or 110V., 1-phase+ neutral + ground, 60Hz. and 3 CV max.

EDIBON Steam Generator (TGV), or steam generator with similar characteristics.

Water supply and drainage.

Computer (PC).

DIMENSIONS & WEIGHTS

TAAC unit :	-Dimensions:1600x570x1500 mm. approx. -Weight : 200 Kg. approx.
Control Interface Box:	-Dimensions: 490x330x310 mm. approx. -Weight: 10 Kg. approx.
PLC Module (PLC-PI):	-Dimensions: 490x330x310 mm. approx. -Weight: 30 Kg. approx.

RECOMMENDED ACCESSORIES

For refilling R134a refrigerant and maintenance, we recommend:

- T/KIT1. Maintenance Kit, containing: vacuum pump, hoses and manometers.
- T/KIT2. Maintenance Kit, containing: leakage detector.
- R134a refrigerant (to be acquired by the customer locally).

AVAILABLE VERSIONS

Offered in this catalogue:

-TAAC. **Computer Controlled Air Conditioning Laboratory Unit.**

Offered in other catalogue:

-TAAB. **Air Conditioning Laboratory Unit.**

* Specifications subject to change without previous notice, due to the convenience of improvements of the product.



C/ Del Agua, 14. Polígono San José de Valderas. 28918 LEGANES. (Madrid). SPAIN.
Phone: 34-91-6199363 FAX: 34-91-6198647
E-mail: edibon@edibon.com WEB site: www.edibon.com

Issue: ED01/09
Date: April/2009

REPRESENTATIVE: