


① Unit: TCLGC. Thermal Conductivity of Liquids and Gases Unit

Always included in the supply:

EDIBON PATENT

Teaching Technique used


SCADA. EDIBON Computer Control System



② Control Interface Box



③ Data Acquisition Board




④ Software for:
- Computer Control
- Data Acquisition
- Data Management



Computer (not included in the supply)

- ⑤ Cables and Accessories
⑥ Manuals



www.edibon.com

Products
Products range
Units
9.-Thermodynamics & Thermotechnics

DESCRIPTION

This unit has been designed to enable students to easily determine the thermal conductivity of liquids and gases.

By the realization of the practices the student can determine the thermal conductivity of any suitable gas or compatible liquid with materials on construction.

The unit consists of an aluminum cylinder that forms the body of the unit which has inside a variable power resistance and temperature sensors, one of them inside the own resistance. It also has a brass jacket that contains the test fluid and the refrigeration water.

The fluid whose thermal conductivity is to be determined, fill the small distance of radial security between the body of heated aluminum and a jacket with cooled water. This space is small enough to avoid the natural convection of the fluid.

The aluminum body is centered in the water jacket thanks to two "o" ring joints, which seal the distance of security, and which allow easily and quick dismantling for cleaning.


The unit incorporates a water regulation valve, which allows the flow of water to be varied; a water flow sensor; and temperature sensors at the water inlet and outlet.

The connection and disconnection of the sensors have been made to be easy and comfortable.

This Computer Controlled Unit, is supplied with the EDIBON Computer Control System (SCADA), including: Control Interface Box + Data Acquisition Board + Computer Control and Data Acquisition Software, for controlling the process and the parameters involved.

PROCESS DIAGRAM AND ELEMENTS ALLOCATION

1 actuator and 7 sensors controlled from any computer, and working simultaneously




OPEN CONTROL
+
MULTICONTROL
+
REAL TIME CONTROL

Note:
ST= Temperature sensor.
SC= Flow sensor.
AVR= Regulation valve.

Items supplied as standard

① TCLGC. Unit:


Bench-top unit.
Anodized aluminium structure and panel in painted steel.
Main metallic elements in stainless steel.
Diagram in the front panel with similar distribution to the elements in the real unit.
Aluminium body (cylinder) with brass jacket that contains the test fluid and the refrigeration water.
Variable heating resistance (in the cylinder), computer controlled, (150 W, temperature max.: 150°C). Resistance power controlled from computer (0-100%). The power is measured by a sensor.
6 Temperature sensors, "J" type.
Water flow sensor.
Water flow regulation valve.
Valves.
Syringe.



TCLGC Unit

② TCLGC/CIB. Control Interface Box :

Control interface box with process diagram in the front panel and with the same distribution that the different elements located in the unit, for an easy understanding by the student.
All sensors, with their respective signals, are properly manipulated from -10V. to +10V. computer output. Sensors connectors in the interface have different pines numbers (from 2 to 16), to avoid connection errors.
Single cable between the control interface box and computer.
The unit control elements are permanently computer controlled, without necessity of changes or connections during the whole process test procedure.
Simultaneously visualization in the computer of all parameters involved in the process.
Calibration of all sensors involved in the process.
Real time curves representation about system responses.
Storage of all the process data and results in a file.
Graphic representation, in real time, of all the process/system responses.
All the actuators' values can be changed at any time from the keyboard allowing the analysis about curves and responses of the whole process.
All the actuators and sensors values and their responses are placed in only one computer screen.
Shield and filtered signals to avoid external interferences.
Real time PID control with flexibility of modifications from the computer keyboard of the PID parameters, at any moment during the process. Real time PID and on/off control for pumps, compressors, resistances, control valves, etc. Real time PID control for parameters involved in the process simultaneously. Proportional control, integral control and derivative control, based on the real PID mathematical formula, by changing the values, at any time, of the three control constants (proportional, integral and derivative constants).
Open control allowing modifications, at any time and in a real time, of parameters involved in the process simultaneously.
Possibility of automatization of the actuators involved in the process.
Three safety levels, one mechanical in the unit, other electronic in control interface and the third one in the control software.



TCLGC/CIB

③ DAB. Data Acquisition Board:


PCI Data acquisition board (National Instruments) to be placed in a computer slot. Bus PCI.
Analog input: Channels= 16 single-ended or 8 differential. Resolution= 16 bits, 1 in 65536.
Sampling rate up to: 250 KS/s (Kilo samples per second). Input range (V)=±10V.
Data transfers=DMA, interrupts, programmed I/O. DMA channels=6.
Analog output: Channels=2. Resolution= 16 bits, 1 in 65536. Max. output rate up to: 833 KS/s.
Output range(V)=±10V. Data transfers=DMA, interrupts, programmed I/O.
Digital Input/Output: Channels=24 inputs/outputs. D0 or DI Sample Clock frequency: 0 to 1 MHz.
Timing: Counter/timers=2. Resolution: Counter/timers: 32 bits.



DAB

④ TCLGC/CCSOF. Computer Control + Data Acquisition + Data Management Software:

Compatible with actual Windows operating systems. Graphic and intuitive simulation of the process in screen.
Compatible with the industry standards.
Registration and visualization of all process variables in an automatic and simultaneously way.
Flexible, open and multicontrol software, developed with actual windows graphic systems, acting simultaneously on all process parameters.
Analog and digital PID control. Menu for PID and set point selection required in the whole work range.
Management, processing, comparison and storage of data.
Sampling velocity up to 250,000 data per second guaranteed. Calibration system for the sensors involved in the process.
It allows the registration of the alarms state and the graphic representation in real time.
Comparative analysis of the obtained data, after the process and modification of the conditions during the process.
Open software, allowing to the teacher to modify texts, instructions.
Teacher's and student's passwords to facilitate the teacher's control on the student, and allowing access at different work levels.
This unit allows that the 30 students of the classroom can visualize simultaneously all results and manipulation of the unit, during the process, by using a projector.



TCLGC/SOF

⑤ Cables and Accessories, for normal operation.

⑥ Manuals: This unit is supplied with 8 manuals: Required Services, Assembly and Installation, Interface and Control Software, Starting-up, Safety, Maintenance, Calibration & Practices Manuals.

***References 1 to 6: TCLGC + TCLGC/CIB + DAB + TCLGC/CCSOF + Cables and Accessories + Manuals are included in the minimum supply, enabling a normal operation.**

Continue...

Complementary items to the standard supply

PLC. Industrial Control using PLC (7 and 8):

⑦ PLC-PI. PLC Module:

Circuit diagram in the front panel.

Front panel:

Digital inputs(X) and Digital outputs (Y) block:

16 Digital inputs, activated by switches and 16 LEDs for confirmation (red).

14 Digital outputs (through SCSI connector) with 14 LEDs for message (green).

Analog inputs block:

16 Analog inputs (-10V. to + 10V.) (through SCSI connector).

Analog outputs block:

4 Analog outputs (-10V. to + 10V.) (through SCSI connector).

Touch screen:

High visibility and multiple functions.

Display of a highly visible status.

Recipe function.

Bar graph function.

Flow display function.

Alarm list.

Multi language function.

True type fonts.

Back panel:

Power supply connector.

Fuse 2A.

RS-232 connector to PC.

Inside:

Power supply outputs: 24 Vdc, 12 Vdc, -12 Vdc, 12 Vdc variable.

Panasonic PLC:

High-speed scan of 0.32 μ sec. for a basic instruction.

Program capacity of 32 Ksteps, with a sufficient comment area.

Free input AC voltage(100 to 240 V AC).

DC input: 16 (24 VDC).

Relay output: 14 (250 VA AC/2 A).

High-speed counter.

Multi-point PID control.

Digital inputs/outputs and analog inputs/outputs Panasonic modules.

Communication RS232 wire, to computer (PC).

⑧ TCLGC/PLC-SOF. PLC Control Software:

For this particular unit, always included with PLC supply.



PLC-PI


Items available on request

⑨ TCLGC/CAL. Computer Aided Learning Software (Results Calculation and Analysis).


⑩ TCLGC/FSS. Faults Simulation System.

Software Main Screens

Main screen



Unit in operation



Note: ST= Temperature sensor. SC= Flow sensor. AR=Heating resistance.

Continue...

Examples of Sensors Calibration screens

NIWriteCal.vi

Analog Input Channel: ST-3

Sensor Name: ST-3

Least Squares Fit

Gain: 107.495 Offset: -6.99457

PTA: 10

Volts: -9.003E- Calibrated: -6.996

ENTER SAVE & EXIT CANCEL

ACTUATORS

AR-1

AA-1

AB-1

0 1 2 3 4 5 6 7 8 9 10

Simultaneous Calibration

Reference Value: 0

Data taken: 0

Reference Select	Sensors	Volts	Calibrated	AT
<input type="checkbox"/>	ST-1	0.2753	28.8346	28.83
<input type="checkbox"/>	ST-2	0.3335	29.7856	29.79
<input checked="" type="checkbox"/>	ST-3	0.331	29.0641	29.06
<input checked="" type="checkbox"/>	ST-4	0.3254	29.5453	29.55
<input checked="" type="checkbox"/>	ST-5	0.3295	29.4276	29.43
<input checked="" type="checkbox"/>	ST-6	0.3458	34.752	34.75
<input type="checkbox"/>		-0.0037	-0.0037	0
<input type="checkbox"/>		-0.004	-0.01826	8.02
<input type="checkbox"/>		3.4768	3.4768	3.48
<input type="checkbox"/>		3.215	291.888	291.89
<input type="checkbox"/>		3.066	3.066	3.07
<input type="checkbox"/>		2.6614	2.6614	2.66
<input type="checkbox"/>		2.4281	2.4281	2.43
<input type="checkbox"/>	SC-1	0.1291	0.1424	0.14
<input type="checkbox"/>	SC-2	0.0104	-0.0211	0.02
<input type="checkbox"/>	AN-1	5.9886	5.9886	5.99

☐ Select all

ENTER DONE

GAIN	OFFSET	r
ST-1	101.705	0.8354
ST-2	98.5001	-3.0594
ST-3	102.291	-4.7913
ST-4	102.262	-3.7268
ST-5	101.438	-3.9967
ST-6	91.5356	3.1025
	1	0
	105.08	-7.5992
	1	0
	92.6831	-6.0846
	1	0
	1	0
	1	0
SC-1	0.784847	0.0411
SC-2	0.9199	-0.0307
AN-1	1	0

EXIT SAVE & EXIT

REQUIRED SERVICES

- Electrical supply: single-phase, 220 V. 50Hz or 110V. 60 Hz.
- Water supply.
- Computer (PC).

DIMENSIONS & WEIGHTS

- | | |
|------------------------|---|
| TCLGC Unit: | -Dimensions: 500 x 400 x 300 mm. approx.
-Weight: 40 Kg. approx. |
| Control Interface Box: | -Dimensions: 490 x 330 x 310 mm. approx.
-Weight: 10 Kg. approx. |
| PLC Module (PLC-PI): | -Dimensions: 490 x 330 x 310 mm. approx.
-Weight: 30 Kg. approx. |

*Specifications subject to change without previous notice, due to the convenience of improvements of the product.



C/ Del Agua, 14. Polígono Industrial San José de Valderas.
28918 LEGANÉS. (Madrid). SPAIN.
Phone: 34-91-6199363 FAX: 34-91-6198647
E-mail: edibon@edibon.com WEB site: www.edibon.com

Issue: ED02/10
Date: May/2010

REPRESENTATIVE:

